




*Crane*  
DENTAL LABORATORY  
Helping You Create Lasting Smiles

or

**tune**<sup>TM</sup>  
clear aligners

powered by G&H Orthodontics



“

*I can truthfully state that the software and aligners compete with the best of them. I've noticed better molar control and torque, more efficient arch expansion, significantly less refinements, faster overall treatment times, and less dependence on attachments.*

**- Dr. David Lach DDS MS (Lach Orthodontic Specialists)\***

\*This statement, provided by Dr. David Lach DDS MS, is solely based on personal, professional experience. This statement is not intended to be a substitute for professional medical advice, diagnosis, or treatment, and does not constitute medical or other professional advice.



# 5%

## Average Refinement Rate

The average refinement rate for tune clear aligners is 5%<sup>1</sup> compared to the < 70% reported by doctors who use other aligner systems<sup>2</sup>




## No Contracts

Pay by the aligner, not by the case, with no case minimums and no annual contracts.

## No Compromise

You maintain complete and direct control of your customer interactions from start to finish



## No Wasted Movements

The SmileStudio algorithm eliminates round-tripping to reduce the number of aligners required and overall treatment time.

## Turnkey Solution

A complete, FDA approved, service from case planning to final case packaging



# UP TO 20%

## Fewer Aligner Stages

Tune clear aligners deliver consistent treatment outcomes in 10-20% fewer stages than other aligner systems<sup>3</sup>

# SmileStudio software

We utilize real-time automatic adjustments and comprehensive tooth-movement algorithms to produce ultimate treatment outcomes.

The intuitive user interface keeps control in the hands of the clinician and offers the flexibility to fully customize every patient's case based on treatment preferences. Attachments, interproximal reduction (IPR), final archform shape, and tooth movement and rotation speed can all be tailored to fit your standards.



Call customer service  
to schedule your  
demo today

## Blending your expertise with the latest technology

Our SmileStudio software was designed with the orthodontist in mind. Doctors can use G&H powered clear aligners to treat a wide variety of malocclusions including:



Crowding



Spacing



Underbite



Overbite



Open Bite



Crossbite


Make tooth position changes and automatically calculate the number of stages needed

Change view angles for individual or both arches

Multiple interactive options help you ensure the ultimate smile

Standard and fully customizable attachments

Make individual tooth adjustments to fine tune your treatment plan




Track changes to aligner stages in real time

The amount, stage and leftover spacing of IPR is automatically calculated

# The latest innovation in aligner plastic technology

% of force

All aligners lose force over time, and the percentage of force delivered drops significantly after initial insertion. Zendura FLX is a patented aligner material designed to deliver best-in-class tooth movement forces over a 7-10 day wear period while also providing exceptional patient comfort and reliable accuracy.\*



0

24

48

72

Time (hours)



Zendura FLX



Leading Aligner Brand A




Leading Aligner Brand B

## Fewer pain points for you, and your patient

By specifically designing the plastic to have less force when initially inserted, but with more elastic recovery, patients have less pain (increased compliance) and more tooth-moving forces over time.\*

# The magic of multiple layers



Zendura FLX features an innovative dual-shell construction that provides exceptional elastic rebound as well as precise fit and strength. The magic is an elastomeric inner layer encased in a rigid outer shell that also grips the teeth firmly and protects the aligner against bruxing or grinding from the opposing dental arch.\*

*Patients all seem to be pleased with the comfort and feel of the plastic as well as the lack of streaking that can be seen with a competitor's trays.*

**- Dr. Zach McCallum, D.D.S\***

*Tried a lot of other brands and Zendura is simply the best. It adapts very well and is very easy to trim and polish.*


**- Dr. Ahmed Mansour, D.D.S\***

\*This information was created by, and is the property of, Zendura Dental. Zendura FLX, Zendura Dental, and all associated claims, data, and trademarks are the property and responsibility of Bay Materials, LLC. The statements made by Dr. McCallum and Dr. Mansour are in reference to the Zendura FLX material and not directly related to any experience with Tune Clear Aligners.



# Pricing you can smile about

**No contracts. No case minimums. No surprises.**



Contracts and case minimums should not dictate your treatment plans – your experience and expertise should. G&H Powered Clear Aligners keep you in complete control – from clinician engagement through delivery of the final case, you deliver exceptional service with the power of the G&H turnkey solution.

**Set-up for success  
(without the fees)**

**Getting started** **G&H POWERED**  
**is easy with** **clear aligners**

\*Introductory pricing is not guaranteed and is subject to change without notice. ARCADLABS inc. clear aligner treatment planning software and manufacturing processes are the property and responsibility of ARCADLABS inc. and distributed in partnership with G&H Orthodontics. Clear aligners are classified as medical devices and must be administered by a licensed medical professional. All aspects of medical treatment, including resulting outcomes of treatment associated with clear aligner trays distributed by G&H Orthodontics are the sole responsibility of the licensed medical professional prescribing treatment.



# EXPRESS

1 - 6 Stages

# \$25\*

*per aligner*

Assuming 2 aligners per stage

# UNLIMITED

7+ Stages

# \$20\*

*per aligner*

Assuming 2 aligners per stage

Private Label  
**PACKAGING**

# \$5\*

*per case*



## The G&H Advantage

- FDA Approved
- Unlimited Refinements
- Full Service Staging Support
- Robust Treatment Planning Software & Support
- Transition & Rotate Simultaneously
- No Annual Case Minimums
- No Setup Costs
- No Case Plan Fees



Made in the USA

[ghorthodontics.com/tune](https://ghorthodontics.com/tune)

1 5% refinement rate based on total refinements reported. Data on file.

2 Based on results from national survey of GPs. Data on file.

3 Based on national average of reported aligner stages. Data on file.

4 Based on national average of reported aligner treatment cost. Data on file.

tune™ clear aligners is a registered trademark of G&H Orthodontics

Available only in select countries

Ver. 1.A.2021